
# Process Safety



# Synthesis of 2,3,6-trimethyl, 4-nitrosophenol


#### Introduction:



2,3,6-trimethylphenol

2,3,6-trimethyl, 4-nitrosophenol

- The synthesis of 2,3,6-trimethyl, 4-nitrosophenol is performed in three steps :
- 1) addition of H2SO4 to 2,3,6-trimethylphenol at 2°C
- 2) addition of sodium nitrite (NaNO2) still at 2°C
- 3) heating up to 25℃ to initiate the reaction



### **Experimental**

The DRC is used with reactors of 250 ml. A stirring at 300 rpm is produced.

The reference and measure reactors are loaded with 138 ml (109g) of ethanol.

**23.03g** of **2,3,6-trimethylphenol** is introduced in the measure reactor.

The solution is stirred for 30 min at  $20^{\circ}$ C and the DRC is cooled down to  $2.1^{\circ}$ C

After this preparation phase, the data collection is started.

The figure hereunder presents the signals corresponding to the whole experiment.



Instrument DRC Evolution (-80 to 150°C).



## www.setaram.com - sales @setaram.com

# Excellence in thermal analysis and calorimetry




A calibration is generally carried out before and after a reaction. The purpose of a calibration is to determine the exchange coefficient UA. UA is mainly a function of filling of the reactor, of the temperature and the nature of the material inside the reactor.

In this example a calibration is carried out between each step. Joule effects are applied in the measure reactor. (5 W during 300 seconds).

The experimental conditions and the calculations are explained for each of the three steps.

## Step 1 - Introduction of 6ml of H2SO4 at 2℃



After the first calibration and when the signal  $\Delta T$  is stable, **6 ml** (10.87g) of sulphuric acid is added with a syringe pump. The whole addition takes about 40 min.

For the calculations, the following values are used:

- Sulfuric acid: Cp = 4.0 J/g.K
- Temperature of the reaction: 2.14℃
- Average temperature of introduction : T<sub>dos</sub> = 22°C
- mass = 109 + 23.03 + 10.87 = 142.9 g
- UA = (UA<sub>1</sub>+ UA<sub>2</sub>)/2 = (3.20+3.03)/2= 3.115 W/℃

(The use of the average value of UA, enables to take into account the volume variation inside the reactor and the corresponding evolution of exchange coefficient)

### Results:

•  $Q_{dos} = -865.9 \text{ J}$ .  $Q_{dos}$  corresponds to the energy necessary to cool the added  $H_2SO_4$  and to bring it from the introduction temperature ( $\approx 22^{\circ}$ ) to the working temperature (2.14 $^{\circ}$ )

• **Q** = 5 889.0 J. It is the total energy of reaction of addition and it is obtained by integration of the peak of reaction and addition of  $Q_{dos}$ .

•  $\Delta T_{adiabatic} = 5889 / (142,9 x Cp) = 41.21 / Cp$  It corresponds to the increase of temperature which would be observed in adiabatic mode during the addition of sulfuric acid.(Cp is the Cp of the mixture ethanol + 2,3,6 trimethylphenol + sulfuric acid.)

## www.setaram.com – sales @setaram.com




## Excellence in thermal analysis and calorimetry

# Process Safety



# Step 2 - Introduction of Sodium Nitrite at 2°C



After the second calibration and when the signal  $\Delta T$  is stable, a solution prepared with 18.03g of Sodium Nitride in 40ml of water is slowly introduced in the measure reactor. The whole addition takes about 1 h 05 min.

For the calculations, the following values are used:


- Sodium nitrite solution: Cp = 4.0 J/g.K
- Temperature of the reaction: 2.14°C
- Average temperature of introduction : T<sub>dos</sub> = 24℃
- UA = UA<sub>2</sub> =3.03W/℃

(In this step, it is considered that the level in the reactor does not change. Consequently, the value of UA obtained aftyer the previous step is used.)

<u>Results:</u>

•  $Q_{dos} = -4549.5J$  Qdos corresponds to the energy necessary to cool the added solution of Sodium Nitrite and to bring it from the introduction temperature ( $\approx 24^{\circ}$ C) to the working temperature (2.14 $^{\circ}$ C)

•  $Q_1 = -21 \ 924.0 \ J$  is the total energy of reaction of addition and it is obtained by integration of the peak of reaction and addition of  $Q_{dos}$ .



## Step 3 - Reaction at 25°C

When the introduction is complete, The temperature of the bath is programmed up to 25°C. *For the calculations, the following value is used:* 

### • UA = 1.69 W/℃

As the heating to  $25^{\circ}$  is rapid, the value of UA determined at the end of the reaction is used for this step.)

### <u>Results</u>

 $Q_2 = -6.687 \text{ J}$  is the total energy of reaction of addition and it is obtained by integration of the peak of reaction at 25°C

## www.setaram.com - sales @setaram.com

