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
Fig.1. HPLC-FIAS-AAS in Taiwan Application

Fig.2. HPLC-ICP/MS in Taiwan Application

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Arsenic Species		LD <sub>50</sub> in Mice (mg/kg)
Arsenious acid	(AsIII)	4.5
Arsenic acid	(AsV)	14-18
Monomethylarsonate	(MMA)	1,800
Dimethylarsinate	(DMA)	2,600
Arsenobetaine	(AsB)	10,000

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

Parameter	Setting
Solvent A	20 mM Ammonium Bicarbonate @ pH 8.5
Solvent B	20 mM Ammonium Sulfate @ pH 7.0
Gradient Profile	6 min at 100% A; Step to 100% B for 12 min
Flow Rate	1.5 mL/min
Column	Anion Exchange, Hamilton PRP-X100, 4.1 mm i.d. x 250 mm, 10 μm
Column Temperature	Ambient
Autosampler Flush Solvent	5% Methanol / 95% DI Type I Water
Sample Injection Volume	100 μL
Re-equilibration Time	15 min
Urine Sample Prep	1:10 with DI Water
Solid Phase Extraction (SPE) Cartridge	3mL Supelclean LC-SAX SPE cartridge (Supelco™) in weak solvent reservoir
Detection	PerkinElmer/SCIEX ELAN 6100
Total Analysis Time	33 min

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Parameter	Setting
Mobile Phase	10 mM Ammonium Nitrate and 10 mM Ammonium Phosphate (dibasic); pH 9.4
Flow Rate	1.5 mL/min
Run Time	10 min
Column	Anion Exchange, Hamilton PRP-X100, 4.1 mm i.d. x 250 mm, 10 μm
Column Temperature	Ambient
Autosampler Flush Solvent	5% Methanol / 95% DI Type I Water
Sample Injection Volume	100 μL
Urine Sample Prep	1:5 with Mobile Phase
Detection	PerkinElmer/SCIEX ELAN DRC II
Total Analysis Time	10 min

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圖四、Isocratic chromatograms of five arsenic species (1000 ng/L each) and a blank acquired with the ELAN DRC monitoring  $AsO^+(m/z=91)$

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